

SOFTWARE FOR PLC

clarity

CLARITY ADVANCED CHROMATOGRAPHY DATA STATION (CDS) WITH SOFTWARE MODULE FOR CONTROL OF ALL SEPARLAB INSTRUMENTS AND DEVICES !

clarity can be used with any PLC system. It perfectly process chromatography data including all quantitative calculations, calibrations and problematic peaks integrating. Together with optional Control modules and Extensions it provides the possibility of complete chromatography handling system. The Control modules provide integrated control of selected instruments, Extensions provide functions for specific separation techniques such as PDA or GPC.

CLARITY DATASTATION SYSTEMS	
Specification	Identification
Clarity Extensible single instrument SW CFR21/GLP ready	C50
Clarity Lite simplified single instrument SW	C40
Specification	Identification
Extension modules for Clarity	
Clarity Add-on single instrument add-on	C55
SST SW module for system suitability test	A22
GPC SW module for GPC/SEC data processing	A28
PDA SW module for PDA data processing	A29
Specification	Identification
Control modules for Clarity	
LC Control SW module for HPLC control	A24
Specification	Identification
Hardware - A/D and D/A converters for Clarity	
Colibrick set external USB, 1 channel, cable	U31
Colibrick set external USB, 2 channels, cable	U32
Colibrick set external USB, 4 channels, cable	U34
Validation Kit peak generator and SW methods for OQ	CVK

